	WEDGLIDY	TEO
100	WEBSURV	IES

Person 7 - View/Update Person Risk Factors - Donated or Received Blood or Tissue

Version #: _	Test Cycle #:	Database: Circle one	DEVO	TEST	PROD
Test Date: _		Role Level: Circle one	USER	QA A	ADMIN
Tester Name:_		Circle one	w/HIV	w/o H	IV
User ID: _		Final Result: Circle on	e PAS	S or F	AIL

The following preconditions must be met:

- 1. Successful logon to Web Surveillance Application
- 2. Search and select party
- 3. Expand History node

Business Rules:

- 1. At any time when the system is preparing to display or saving information, the system may show a "Building ... page..." or "Please wait..." message in the right section.
- 2. Comments/Donated / Received cannot be edited. They can only be removed prior to being saved. Punctuation and upper/lower case is allowed for Donated / Received.
- 3. Users may add and edit information before it is added to the database, but cannot update after the person records have been saved. Some items may only be removed.
- 4. QA and Admin users may add and edit information before and after the person records have been saved. Some items may only be removed.
- 5. Error messages are shown in red above the section of the screen.
- 6. At any time when you click Cancel at bottom of screen, the screen is reloaded.
- 7. Donating Site can only be entered when Donated.
- 8. Receiving Organization can only be entered when Donated.
- 9. Healthcare Facility can only be entered when Received.

Procedure #1	Expected System Response		Fail
Click on the Donated / Received Node	The Donated / Received Screen is		
	returned.		
Procedure #2	Expected System Response	Pass	Fail
Select "D" in Donated/Received.	The screen is expanded to include		
	donating site and receiving organization.		
Procedure #3	Expected System Response	Pass	Fail
Enter date, blood or tissue type, and	The Search Organizations screen is		
donating site. Click search for Receiving	returned.		
Organization.			
	T	-	
Procedure #4	Expected System Response	Pass	Fail
Enter Organization name or select	The Search Organizations screen is		
organization type and click Search	returned with results.		
Procedure #5	Evacated System Bearense	Pass	Fail
	Expected System Response	F a 5 5	ган
Select an organization by clicking on the	The screen is returned with the		
name link.	receiving organization shown.		

◯ WEBSURV TEST

Procedure #6	Expected System Response	Pass	F ail
Click Add to List	Expected System Response Validation failed, messages	rass	raii
WW Click Add to List	returned.		
	Validation passed, the record is		
	shown in the grid.		
	Shown in the gha.	1	
Procedure #7	Expected System Response	Pass	Fail
Select "R" in Donated/Received	The screen is expanded to include		
	Health Care Facility.		
Procedure #8	Expected System Response	Pass	F ail
Enter date, blood or tissue type, and	The Search Organizations screen is		
Click search for Healthcare Facility.	ty. returned.		
Procedure #9	Evacated System Response	Pass	F ail
Senter organization name or select	Expected System Response	rass	rall
	The Search Organization screen is returned with results.		
organization type and click search.	returned with results.		
Procedure #10	Expected System Response	Pass	Fail
Select an organization by clicking on the	The Healthcare Facility is shown on	газз	ran
name link.	the screen.		
name iink.	the screen.	<u> </u>	
Procedure #11	Expected System Response	Pass	Fail
Click Add to List	Validation failed, messages	1 400	1 (
o oner rad to blet	returned.		
	Validation passed, the record is		
	added to the grid.		
	9		I .
Procedure #12	Expected System Response	Pass	F ail
Click Edit on a record.	The information is shown in the		
	entry fields.		
Procedure #13	Evented System Response	Doos	Fail
	Expected System Response	Pass	Fail
Senter/change information and click	The changed information is shown in		
Replace.	the grid.		
Procedure #14	Expected System Response	Pass	Fail
Click Remove next to record.	A confirm remove message is	1 433	• an
Silek Remove flext to record.	returned.		
L			1
Procedure #15	Expected System Response	Pass	Fail
♦ Click OK	The information is removed from		
	grid.		
Procedure #16	Expected System Response	Pass	Fail
Enter comments and click Add to List	The comment is added to the grid.		
Procedure #17	Expected System Response	Pass	Fail
V Click Save	The View Donated / Received	. 433	. un
J Chief Cave	screen is returned.		
	Jordon is returned.		

Test is Complete.

Print Date: 9/24/2008